



Quality Control Ranges

**MILLIPLEX[®] MAP Mouse Amyloid Beta
Magnetic Bead Panel
Kit Cat# MABMAG-83K
Cat# MABMG-6083**

<u>Analyte</u>	<u>QC Level</u>	<u>Lot #</u>	<u>Expected Range</u>	<u>Units</u>
A β 40	Control 1	MAB83-104	37 - 78	pg/mL
	Control 2	MAB83-204	262 - 544	pg/mL
A β 42	Control 1	MAB83-104	70 - 146	pg/mL
	Control 2	MAB83-204	487 - 1011	pg/mL

MABMG-6083: MAB83-104.204 (06 August 2024)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA