



QUALITY CONTROL RANGES

HUMAN GHRELIN (TOTAL) ELISA CATALOG # EZGRT-89K

STANDARD, LYOPHILIZED CATALOG # E8089-K

STANDARD	HUMAN GHRELIN VALUE
LOT NUMBER: EGT-13K	4348 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **2.0mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6089-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	TGR-113	243 – 504	pg/mL
QC 2	TGR-213	1047 – 2174	pg/mL

E6089-K: TGR-113.213 (17 July 2020)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA.