



RABBIT ANTI-HUMAN CYTOCHROME P450 ENZYME CYP1A2 POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB1248	QUANTITY:	100 µL
LOT NUMBER:			
SPECIFICITY:	Reacts with human cytochrome P450 enzyme CYP1A2 in hepatic microsomal fraction. No cross-reactivity with other cytochrome P450 enzymes. By Western blot the antibody recognizes a protein with a molecular weight of approximately 54 kDa.		
IMMUNOGEN:	Synthetic peptide		
APPLICATIONS:	Western blot: 1:1,000-1:4,000 Immunohistochemistry on frozen and formaldehyde treated sections. Optimal working dilutions must be determined by end user.		
SPECIES REACTIVITY:	Human and rat.		
FORMAT:	Rabbit antiserum		
PRESENTATION:	Liquid. No preservatives.		
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 12 months. Avoid repeated freeze/thaw cycles.		
RELATED REFERENCES:	<i>Federation Proc.</i> (1987) 46 :2563-2566. <i>Adv. Enzymol.</i> (1987) 60 :35-87. <i>DNA and Cell Biology</i> (1993) 12 :1-51. <i>Biochemical Pharmacology</i> (1995) 49 :39-47.		

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 - 2011: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permission in writing.