



RABBIT ANTI- KISSPEPTIN 10 (kp10) POLYCLONAL ANTIBODY

CATALOG NUMBER:	AB9754
LOT NUMBER:	
QUANTITY:	100 µL
SPECIFICITY:	Kisspeptin 10 (kp10, metastin). No cross reactivity with GnRH, Galatin, SP, NPY, CRH, alphaMSH, SRIF or PrRP by RIA.
IMMUNOGEN:	Peptide from mouse kisspeptin 10. The sequence is common to all kisspeptin isoforms.
APPLICATIONS:	Immunohistochemistry on 4% paraformaldehyde fixed tissue: 1:500-1:1,000 (PAP and ABC with NiDAB detection system). RIA: 1:500-1:10,000 Optimal working dilutions must be determined by the end user.
SPECIES REACTIVITY:	Rat, mouse and ovine. The antibody showed a very low level of reactivity with human kp10. Other species have not been tested.
FORMAT:	Rabbit serum.
PRESENTATION:	Liquid containing 50% glycerol, Ph 7.4, Phosphate buffer 0.01M.
STORAGE/HANDLING:	Maintain at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.
REFERENCE:	Franceschini, I., et al., <i>Neuroscience Letters</i> (2006) 401 :225-230.

Important Note: During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.

FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.