

3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

# **ProductInformation**

## Minocycline hydrochloride

Product Number **M 9511** Storage Temperature 2-8 °C

## **Product Description**

Molecular Formula: C<sub>23</sub>H<sub>27</sub>N<sub>3</sub>O<sub>7</sub> • HCl

Molecular Weight: 493.9 CAS Number: 13614-98-7

Synonyms:  $[4S-(4\alpha,4a\alpha5,5a\alpha,12a\alpha)]-4,7$ -bis(dimethylamino)-14,4a,5,5a,6,11,12a-octahydro-3,10,12,12a-tetrahydroxy-1,11-dioxo-2-naphthacenecarboxamide; 7-dimethylamino-

6-demethyl-6-deoxytetracycline<sup>1</sup>

Minocycline is a semi-synthetic tetracycline derivative that has a spectrum of antibacterial activity similar to tetracycline. It is active against streptococci, enterobacteria, and some mycobacteria, and also against such species as *Staphylococcus aureus*, *Neisseria meningitidis*, *Acinetobacter*, *Bacteroides*, *Haemophilus*, and *Nocardia*. Minimum inhibitory concentrations (MIC) can range from 0.06 - 1 μg/ml for the most sensitive organisms, and from 4 - 12.5 μg/ml for moderately sensitive organisms.<sup>1,2</sup> A detailed review of minocycline has been published.<sup>3</sup>

Minocycline has been shown to inhibit tumor growth in transgenic mice, in combination with AGM-1470 and interferon  $\alpha/\beta.^4$  Minocycline has also been demonstrated to inhibit angiogenesis in rabbit cornea in the presence of the VX2 carcinoma.  $^5$  The inhibitory activity of minocycline against several metalloproteinases has been studied.  $^6$ 

An HPLC assay for the analysis of minocycline in serum and tissue has been published.<sup>7</sup>

#### **Preparation Instructions**

This product is soluble water (50 mg/ml), with heat as needed, yielding a clear, yellow to amber solution. It is also soluble in DMSO (7 mg/ml).

## Storage/Stability

Stock solutions of this product (1 mg/ml, 0.1 N HCl) may be stored for two days at 4 °C.

#### References

- 1. The Merck Index, 12th ed., Entry# 1188.
- 2. Martindale The Extra Pharmacopoeia, 31st ed., Reynolds, J. E. F., ed., Royal Pharmaceutical Society (London, UK: 1996), pp. 829-830.
- 3. Zbinovsky, V., & Chrekian, G. P., in Anal. Prof of Drug Substances, **6**, Academic Press (New York, NY: 1977), pp. 323-339.
- 4. Parangi, S., et al., Antiangiogenic therapy of transgenic mice impairs *de novo* tumor growth. Proc. Natl. Acad. Sci. USA, **93**, 2002-2007 (1996).
- Tamargo, R., et al., Angiogenesis inhibition by minocycline. Cancer Res., 51, 672-675 (1991).
- Gilbertson-Beadling, S., et al., The tetracycline analogs minocycline and doxycycline inhibit angiogenesis in vitro by a non-metalloproteinasedependent mechanism. Cancer Chemother. Pharmacol., 36(5), 418-424, and 37(1-2), 194 (erratum) (1995).
- 7. Wrightson, W. R., et al., Analysis of minocycline by high-performance liquid chromatography in tissue and serum. J. Chromatogr. B Biomed. Sci. Appl., **706(2)**, 358-361 (1998).

GCY/NSB 12/03

### **Precautions and Disclaimer**

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.

For Laboratory Use Only. Not for drug, household or other uses.