



**QUALITY CONTROL RANGES**  
**MILLIPLEX® MAP Multi-Species Hormone**  
**Kit Catalog# MSHMAG-21K**

**Testosterone Standard Catalog # MSHTES-8180**  
**Lot# MSHTES-1K**

<b>Hormone</b>	<b>Standard #7 Value Assignment</b>
Testosterone	150 ng/mL

**Control Catalog#**  
**MSHTES-6180-1 and MSHTES-6180-2**

<b><u>Hormone</u></b>	<b><u>QC Level</u></b>	<b><u>Lot #</u></b>	<b><u>Expected Range</u></b>	<b><u>Units</u></b>
Testosterone	Control 1	MSHTES-101	1.5 – 3.1	ng/mL
	Control 2	MSHTES-201	17 – 35	ng/mL

MSHTES-6180: MSHTES-101.201 (28 April 2025)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • [emdmillipore.com](http://emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA