

80331 EC 0157:H7 ChromoSelect Selective Supplement I

Recommended for isolation and easy detection of Escherichia coli 0157:H7 from food samples.

Composition

(per vial, sufficient for 500 ml medium)

Distilled water 5 ml

Novobiocin 15 mg

Potassium tellurite 1.5 mg

Directions

Aseptically add the contents of 1 vial to 495 ml of sterile, molten Enrichment ChromoSelect Broth (Base) for EC 0157:H7 (Cat. No. 80330), cooled to 45-50 °C. Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

