

Product Information

Anti-MTMR2
antibody produced in chicken,
affinity isolated antibody

Catalog Number **GW21545**
Formerly listed as GenWay Catalog Number 15-288-21545,
Myotubularin-related protein 2 Antibody.

Storage Temperature Store at -20 °C

Synonyms: Myotubularin-related protein 2 isoform 1, EC 3.1.3.-

Product Description

Interacts with This protein is a member of the myotubularin family and encodes a putative tyrosine phosphatase. Mutations in this gene are a cause of Charcot-Marie-Tooth disease type 4B. an autosomal recessive demyelinating neuropathy. Multiple alternatively spliced transcript variants have been found. but the biological validity of some variants has not been determined.

NCBI Accession number: NP_057240.3
Swiss Prot Accession number: Q13614

Gene Information: Human .. MTMR2 (8898)

Immunogen: Recombinant protein Myotubularin-related protein 2 isoform 1

Immunogen Sequence: GI # 44680154, sequence 531 - 591

The product is a clear, colorless solution in phosphate buffered saline, pH 7.2, containing 0.02% sodium azide.

Species Reactivity: Human

Tested Applications: WB

Recommended Dilutions: Recommended starting dilution for Western blot analysis is 1:500, for tissue or cell staining 1:200.

Note: Optimal concentrations and conditions for each application should be determined by the user.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Due to the sodium azide content a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one week. For extended storage, store in -20 °C freezer in working aliquots. Repeated freezing and thawing, or storage in "frostfree" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

LPG,MAM12/09-1