



IHC SELECT® BIOTINYLATED SECONDARY ANTIBODY GOAT ANTI-RAT IgG

Catalog Number:	21543
Volume:	15 mL
Format:	Goat Anti-Rat IgG in PBS containing carrier protein, with 0.2% Tween 20 and 0.1% ProClin 300 as preservative.
Application:	<p>Prediluted for use in Immunohistochemical staining procedures. For example:</p> <ul style="list-style-type: none">○ Blocking Reagent○ Primary Antibody (Rat IgG) or Negative Control Reagent○ Biotinylated Secondary Antibody (Gt X Rat IgG)○ Streptavidin-Enzyme Conjugate○ Chromogen○ Counterstain <p>For best results use this reagent in conjunction with other IHC Select™ detection reagents. A complete list of reagents is available at www.millipore.com.</p>
Storage / Handling:	When stored at 2-8°C, this reagent is stable up to the expiration date printed on the label. Do not freeze or expose to elevated temperatures.

General Purpose Reagent