

Product Information

68835 1-Methyl-piperidine Solution **1-Methyl-piperidine-Lösung**

Dear customer

To adapt our bottle to the most popular ABI protein sequencer, please use the appropriate adapter (product number 17853). This adapter is reusable and will be sent out free of charge following your order of 1-methyl-piperidine solution (product number 68835).

If you have not received the adapter up till now, please order it from our local subsidiary.

Sehr geehrter Kunde

Um die gelieferte Flasche an die weit verbreiteten ABI-Sequenatoren anzuschliessen, benötigen Sie einen Adapter (Produkt-Nr. 17853). Dieser Adapter ist wiederverwendbar und wird Ihnen kostenlos zur Verfügung gestellt.

Falls Sie noch keinen Adapter erhalten haben, fordern Sie ihn bitte bei unseren Niederlassungen an.

Precautions and Disclaimer:

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.