

14121 Lactic Acid Supplement, modified

A selective supplement recommended for the isolation of lactic acid bacteria in beer and brewing process. It stimulates lactic acid bacteria and inhibits Gram-negative organisms, yeasts and moulds.

Composition

(per vial, sufficient for 500 ml medium):

Sorbitan mono-oleate	5.0 ml
2-phenyl ethanol	1.5 g
Amphotericin	2.5 mg

Directions

Aseptically add the content of 1 vial to 500 ml sterile, molten Rakay-Ray Agar (Cat. No. 02538). Mix well and pour into sterile petriplates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

