

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Methanol

Product Number M 1770 Store at Room Temperature

Product Description

Molecular Formula: CH₃OH Molecular Weight: 32.04 CAS Number: 67-56-1 Density: 0.7915 g/ml (20 °C)¹ Boiling point: 64.7 °C (760 torr)¹

Synonyms: methyl alcohol, methyl carbinol¹

This product is as suitable for use as a protein sequencing reagent.

Methanol is a widely used polar organic solvent in chemistry, biochemistry, and molecular biology. It is a common solvent in HPLC and mass spectrometry (MS) of small molecules and of biomolecules such as oligonucleotides and proteins. ^{2,3,4,5} Methanol is also utilized in high performance capillary electrophoresis.6 Protocols which utilize methanol have been published on the isolation of hydrophobic and membrane proteins from macrophages, plants, and bile. 7,8,9 The use of methanol in studies of sequencing gels has been described.¹⁰

Industrial applications of methanol include the manufacture of formaldehyde and methyl esters of organic and inorganic acids. Methanol is used as a solvent and solvent adjuvant for polymers, and as a softening agent for pyroxylin plastics.

Precautions and Disclaimer

For Laboratory Use Only. Not for drug, household or other uses.

Preparation Instructions

This product is miscible with water, ethanol, ether, benzene, and most other organic solvents.

References

- 1. The Merck Index, 12th ed., Entry# 6023.
- 2. Cech, N. B., and Enke, C. G., Practical implications of some recent studies in electrospray ionization fundamentals. Mass Spectrom. Rev., **20(6)**, 362-387 (2001).

- 3. Fountain, K. J., et al., Analysis of native and chemically modified oligonucleotides by tandem ion-pair reversed-phase high-performance liquid chromatography/electrospray ionization mass spectrometry. Rapid Commun. Mass Spectrom., **17(7)**, 646-653 (2003).
- High-Performance Liquid Chromatography of Peptides and Proteins: Separation, Analysis, and Conformation. Mant, C. T., and Hodges, R. S., eds., CRC Press (Boca Raton, FL: 1991), pp. 292, 310, 442, 644.
- 5. Piette, V., et al., Enantiomer separation of Nprotected amino acids by non-aqueous capillary electrophoresis and high-performance liquid chromatography with tert-butyl carbamoylated quinine in either the background electrolyte or the stationary phase. J. Chromatogr. A., 987(1-2), 421-427 (2003).
- 6. High Performance Capillary Electrophoresis: Theory, Techniques, and Applications. Khaledi, M. G., ed., John Wiley & Sons (New York, NY: 1998), pp. 527, 529-531, 536, 539, 541, 546-550,
- 7. Simoes-Barbosa, A., et al., Solubilization of delipidated macrophage membrane proteins for analysis by two-dimensional electrophoresis. Electrophoresis, 21(3), 641-644 (2000).
- Seigneurin-Berny, D., et al., Technical Advance: Differential extraction of hydrophobic proteins from chloroplast envelope membranes: a subcellularspecific proteomic approach to identify rare intrinsic membrane proteins. Plant J., 19(2), 217-228 (1999).
- 9. Stark, M., et al., Isolation and characterization of hydrophobic polypeptides in human bile. Eur. J. Biochem., 266(1), 209-214 (1999).
- 10. Molecular Cloning: A Laboratory Manual, 3rd ed., Sambrook, J., and Russell, D. W., CSHL Press (Cold Spring Harbor, NY: 2001), pp. 12.90-12.92, A1.20.

GCY/RXR 4/03