

18038 Fraser Selective Supplement

Recommended by ISO Committee for *Listeria monocytogens* from food, animal feed etc.

Composition:

(per vial, sufficient for 500 ml Fraser Broth or 1000 ml Half Fraser Broth)

Nalidixic acid	10.0 mg
Acriflavine Hydrochloride	12.5 mg

Directions:

Rehydrate the contents of 1 vial aseptically with 10 ml of sterile distilled water, mix well. Aseptically add 1 vial to 1000 ml sterile Fraser Broth Base (Cat. No. 69198) for primary enrichment (Half Fraser Broth) or to 500 ml sterile, Fraser Broth Base (Cat. No. 69198) for secondary enrichment (Fraser Broth). Mix well and dispense as desired.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

