



QUALITY CONTROL RANGES

Human C-Peptide ELISA

Catalog #: EZHCP-20K

Control Catalog #: E6020-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	HCP-122	0.45 - 0.94	ng/mL
QC 2	HCP-222	3.0 – 6.3	ng/mL

E6020-K: HCP-122.222 (09 December 2025)

SigmaAldrich • [Sigmaaldrich.com/support](https://sigmaaldrich.com/support)