



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sia.com
sigma-aldrich.com

Product Information

Anti-Potassium Channel K_{IR}3.3 (GIRK3)

Developed in Rabbit, Affinity Isolated Antibody

Product Number **P 8247**

Product Description

Anti-Potassium Channel K_{IR}3.3 (GIRK3) was developed in rabbit using a GST fusion protein with sequence (C)RLDA HLYWS IPSRL DEKV, corresponding to residues 344-361 of rat K_{IR}3.3 (GIRK3) as the immunogen. This sequence is identical in mouse and human. The antibody was affinity isolated on immobilized immunogen.

Anti-Potassium Channel K_{IR}3.3 (GIRK3) recognizes K_{IR}3.3 (GIRK3) by Western blotting of rat brain membranes.

The action of potassium (K⁺) channels is regulated by voltage, calcium and a variety of neurotransmitters. Each subfamily generally consists of a primary pore forming α subunit that is associated with several regulatory subunits.¹ To date, some 70 different genes that encode the α subunits of K⁺ channels have been identified. Recently, the crystal structure of the K⁺ channels has been identified.²

The vast family of K⁺ channels has been subdivided into the three main subfamilies: the 2 TM, 4 TM and 6 TM K⁺ channels.³ The G-protein-activated inwardly rectifying potassium channels (GIRKS) are members of the 2 TM family, also known as inwardly-rectifying potassium (K_{IR}) channels. Inward rectifiers have two main physiological roles: to mediate transport across the cell membrane and to stabilize the resting membrane potential near the potassium equilibrium potential.⁴ Four GIRKs, referred to as K_{IR}3.1-3.4, have been identified in mammals.⁵

Reagent

The antibody is supplied as lyophilized powder from phosphate buffered saline, pH 7.4, containing 1% BSA 5% sucrose, and 0.025% sodium azide as preservative.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazards and safe handling.

Preparation Instructions

Reconstitute with 0.05 ml or 0.2 ml deionized water, depending on package size. Further dilutions should be made using a carrier protein such as BSA (1%).

Storage/Stability

Lyophilized powder can be stored at room temperature for several weeks. For extended storage, it should be stored at -20 °C. The reconstituted solution can be stored at 2-8 °C for up to 2 weeks. For longer storage, freeze in working aliquots. Avoid repeated freezing and thawing, and storage in "frost-free" freezers. Centrifuge before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

The recommended working dilution is 1:200 for immunoblotting.

Note: In order to obtain best results in different techniques and preparations we recommend determining optimal working concentration by titration test.

References

- Alexander, S.P., et al., Br. J. Pharmacol., **141**, Suppl 1:S1-S126 (2004).
- MacKinnon, R., FEBS Letters, **555**, 62-65 (2003).
- Gutman, G.A., et al., Pharmacol. Rev., **55**, 583-586 (2003).
- Yamada, M. et al., Pharmacol. Rev., **50**, 723-760 (1998).
- Mark, M.D. and Herlitze, S., Eur. J. Biochem., **267**, 5830-5836 (2000).

MCT/PHC 11/04

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.