



QUALITY CONTROL RANGES

Rat/Mouse C-Peptide 2 ELISA
Catalog #: EZRMCP2-21K

QC Catalog #: E6021-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RCP-114	53 – 110	pM
QC 2	RCP-214	494 – 1026	pM

E6021-K: RCP-114.214 (04 January 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA



E6021-K: RCP-114.214 (04 January 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA