

3050 Spruce Street Saint Louis, Missouri 63103 USA Telephone 800-325-5832 • (314) 771-5765 Fax (314) 286-7828 email: techserv@sial.com sigma-aldrich.com

ProductInformation

Anti-phospho Tau (pSer²⁰²)

Developed in Rabbit, Affinity Isolated Antibody

Product Number T 8069

Product Description

Anti-phospho-Tau (pSer²⁰²) is developed in rabbit using a synthetic phosphopeptide derived from the region of tau that contains serine 202 as immunogen. The serum is affinity purified using epitope-specific affinity chromatography. The antibody is preadsorbed to remove any reactivity towards a non-phosphorylated tau. Anti-phospho-Tau (pSer²⁰²) recognizes human, mouse and rat tau (pSer²⁰²) (45-68 kDa). It has been used in immunoblotting and immunohistochemistry applications.

Tau is a microtubule-associated phosphoprotein (MAP), localized in neuronal axons. It promotes tubulin polymerization and stabilizes microtubules.¹ The biological activity of tau is regulated by its degree of phosphorylation.^{1,2} Hyperphosphorylated tau is the major protein of the paired helical filaments (PHFs), which make up the pathological neurofibrillary tangles of Alzheimer's disease (AD). The PHFs are also found in the lesions of other central nervous system disorders.^{3,4}

Tau phosphorylation involves numerous kinases: glycogen synthase kinase 3β (GSK- 3β), MARK kinase, MAP kinase, protein kinase A and C, cyclin-dependent kinase 5 (Cdk5), p38 kinase, c-Jun N-terminal kinase, and casein kinase II.^{1,2,5,6,7} Combined tau protein kinase II (TPKII), which consists of Cdk5 and GSK- 3β , is the most potent phosphorylation agent indirectly involved in the regulation of the phosphorylation state of tau in neuronal cells.^{6,8} In addition, tau is phosphorylated *in vitro* by *o*smotic cellular stress, which activates the stress-activated protein kinases (SAPKs).

To date, a total of 25 abnormal phosphorylation sites have been identified on hyperphosphorylated tau in AD brain.¹⁰ Normal tau has approximately eight phosphorylation sites. The abnormal phosphorylation occurs usually on serine and threonine residues. Specifically, TPKII phosphorylates serines 202 and 404. GSK-3 β transfection phosphorylates serines 199, 202, 235, 396, 404 and 413, and threonines 205 and 231. These sites are among the major abnormal phosphorylation sites of tau.¹¹ Phosphorylation on these sites reduce the ability of given tau species to promote microtubule self-assembly.^{11,12}. Okadaic acid increases phosphorylation at threonine 231 and serines 235, 396 and 404. Phosphorylated serine 422 was found in the biopsies of brains from patients with Down syndrome, amyotropic lateral sclerosis, corticobasal degeneration, and Pick's disease. It was absent from control group of normal brains.13

The opposite process, tau dephosphorylation, is controlled by different protein phosphatases expressed in neurons. Protein phosphatases PP2A and PP2B efficiently dephosphorylate tau *in vitro* and restore biological activity in the assembly of microtubules.^{3,10,14}

Recently it was discovered that propyl isomerase (Pin 1) interacts with tau hyperphosphorylated on threonine 231 and restores the ability of tau to bind to microtubules.

Reagent

Anti-Phospho-Tau (pSer²⁰²) is supplied as a solution in Dulbecco's phosphate buffered saline (without Mg²⁺ and Ca²⁺), pH 7.3, with 50% glycerol, 1.0 mg/ml BSA (IgG and protease free) and 0.05% sodium azide.

Precautions and Disclaimer

Due to the sodium azide content, a material safety data sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution. Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

Store at -20 °C. Due to the presence of 50% glycerol the antibody will remain in solution. For extended storage, centrifuge the vial briefly before opening and prepare working aliquots. The antibody is stable for at least six months when stored appropriately. Working dilutions should be discarded if not used within 12 hours.

Product Profile

A recommended working dilution of 1:1000 is determined by immunoblotting using cell extracts from African green monkey kidney (CV-1). Data show the specificity of Anti-Phospho-Tau (pSer²⁰²) for the phosphorylated protein since no blocking is observed in the presence of the non-phosphorylated immunizing peptide.

Note: In order to obtain best results in different techniques and preparations we recommend determining optimal working concentration by titration test.

References

- 1. Haque, N., et al., Brain Res., 838, 69-77 (1999).
- 2. Jenkins, S.M., et al., Biochem. J., **345**, 263-270 (2000).
- 3. Gong, C.X., et al., J. Biol. Chem., **275**, 5535-5544 (2000).
- 4. Goedert, M., Nature, **399**, 739-740 (1999).
- 5. Maccioni, R.B., et al., Eur. J. Biochem., **268**, 1518-1527(2001).
- Arioka, M. et al., J. Neurochem., 60, 461-468 (1993).
- Anderton, B.H., et al., Biochem. Soc. Symp., 67, 73-80 (2001).
- 8. Michel, G., et al., Biochem. Biophys. Acta, **1380**, 177-182 (1998).
- 9. Lundt, E.T., et al., J. Neurochem., **76**, 1221-1232 (2001).
- 10. Gong, C.X., et al., J. Neurochem., **62**, 803-806 (1994).
- 11. Takashima, A., et al., Neurosci. Res., **31**, 317-323 (1998).
- 12. Utton, M.A., et al., Biochem. J., **323**, 741-747 (1997).
- 13. Bussiere, T., et al., Acta Neuropathol. (Berlin), **97**, 221-230 (1999).
- 14. Wang, J.Z., et al., Mol. Brain Res., **38**, 200-208 (1996).

AH/JK 7/2004

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.