



QUALITY CONTROL RANGES

MILLIPLEX® MAP Mouse CVD Magnetic Bead Panel 1
Kit Catalog # MCVD1MAG-77K
Control Catalog # MCVDMG1-6077-1
Lot# MCVD1-109 and MCVD1-209

<u>Analyte</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
E-Selectin	Control I	1.3 - 2.6	ng/mL
	Control II	12 - 26	ng/mL
ICAM-1	Control I	0.18 - 0.37	ng/mL
	Control II	1.9 - 4.0	ng/mL
Pecam-1	Control I	0.55 - 1.1	ng/mL
	Control II	5.8 - 12	ng/mL
P-Selectin	Control I	11 - 22	ng/mL
	Control II	106 - 220	ng/mL
PAI-1Total	Control I	0.59 - 1.2	ng/mL
	Control II	5.5 - 11	ng/mL
proMMP-9	Control I	0.30 - 0.63	ng/mL
	Control II	3.3 - 6.8	ng/mL
Thrombomodulin	Control I	0.53 - 1.1	ng/mL
	Control II	6.1 - 13	ng/mL

MCVDMG1-6077-1: MCVD1-109.209 (15 MAY 2025)

[SigmaAldrich.com/support](https://www.sigmaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.