

Supelco

www.sigmaaldrich.com

1.14403

MQuant®

Iron Test

Eisen-Test

Test Fer

Test Hierro

Color card

Farbkarte

Carte colorimétrique


Tarjeta colorimétrica

Fe




© 2024

1.



2.



7.69144.0304-XXXXXX 03/2024

Merck KGaA, 64271 Darmstadt, Germany,
Tel. +49(0)6151 72-2440
www.sigmaaldrich.com/colorimetry

EMD Millipore Corporation, 400 Summit Drive
Burlington MA 01803, USA, Tel. +1-978-715-4321

MilliporeSigma Canada Ltd., 2149 Winston Park Dr,
Oakville, Ontario, L6H 6J8, Canada,
Phone: +1 800-565-1400

0 0.01 0.02 0.03 0.04 0.06 0.08 0.10 0.15 0.20 mg/l Fe



20 ml



5 x Fe-1



3 min

