



QUALITY CONTROL RANGES

MOUSE ADIPONECTIN ELISA CATALOG #: EZMADP-60K STANDARD, LYOPHILIZED CATALOG #: E8060-K

STANDARD LOT NUMBER: EMA60-23K	MOUSE ADIPONECTIN VALUE 153 ng/mL
-----------------------------------	---

The exact value of lyophilized standard is shown above.

Reconstitute with **1 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #: E6300-K

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	QC-123	2.8 – 5.9	ng/mL
QC 2	QC-223	10 – 22	ng/mL

E6300-K: QC-123.223 (10 September 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.

The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.