

## Product Information

## Interleukin-2 Human

Recombinant, Expressed in *Pichia pastoris*, Suitable for Cell culture.

**I7908**

### Description

Interleukin-2 Human, Recombinant, is provided as a DNA sequence encoding the mature human IL-2 protein sequence with cysteine 125 mutated to alanine and expressed in the yeast, *Pichia pastoris*. This mature human IL-2 contains 133 amino acid residues and has a predicted molecular mass of ~15 kDa.

Interleukin-2 (IL-2), also known as T Cell Growth Factor (TCGF), is an immunomodulatory factor produced by certain subsets of T lymphocytes.<sup>1</sup> This lymphokine promotes long term growth of activated T cells and related cell types. Interleukin-2 plays a role in the activation and proliferation of NK cells, induces  $\gamma$ -interferon and B cell growth factor secretion,<sup>2-5</sup> and modulates the expression of the IL-2 receptor.<sup>6</sup> Interleukin-2 has been isolated from various cell types<sup>7-9</sup> and produced by recombinant DNA technology.<sup>10</sup>

### Reagent

Supplied as a lyophilized powder from a 0.2  $\mu$ m filtered solution in 10 mM Phosphate Buffer, pH 7.0, containing 0.3% Human Serum Albumin.

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

Store at -20 °C. After reconstitution, store working aliquots at -20 °C. Repeated freezing and thawing is not recommended and will result in decreased biological activity.

### Preparation Instructions

Reconstitute the contents of the vial in 50 mM Phosphate Buffer, pH 7.0, containing 0.1% Bovine Serum Albumin or Human Serum Albumin (0.1 mL). Upon reconstitution, the cytokine can be stored at 2-8 °C for up to one month without detectable loss of activity. For prolonged use, aliquot and store at -20 °C.

### Product Profile

The biological activity of Recombinant Human Interleukin-2 is determined in a proliferation assay using an IL-2 dependent murine cytotoxic T cell line, CTLL-2.<sup>11</sup>

The ED<sub>50</sub> is defined as the effective concentration of growth factor that elicits a 50% increase in cell growth in a cell-based bioassay.

Purity:  $\geq 98\%$  (SDS-PAGE and HPLC).

Endotoxin: <10 EU (endotoxin units)/ $\mu$ g of cytokine (LAL method).

## References

1. Smith, K. A., Science, 240, 1169 (1988).
2. Morgan, D. A., et al., Science, 193, 1007 (1976).
3. Ortaldo, J. R., et al., J. Immunol., 133, 779 (1984).
4. Farrar, J. J., et al., Immunol. Rev, 63, 129 (1982).
5. Inaba, K., et al., J. Exp. Med., 158, 2040 (1983).
6. Smith, K. A. and Cantrell, D. A. Proc. Natl. Acad. Sci. USA, 82, 864 (1985).
7. Welte, K., et al., J. Exp. Med., 156, 454 (1982).
8. Kniep, E. M., et al., Eur. J. Biochem., 143, 199 (1984).
9. Robb, R. J., et al., Proc. Natl. Acad. Sci., USA, 80, 5990 (1983).
10. Taniguchi, T., et al., Nature, 302, 305 (1983).
11. Gillis, S., et al., J. Immunol., 120, 2027 (1978).

## Notice

We provide information and advice to our customers on application technologies and regulatory matters to the best of our knowledge and ability, but without obligation or liability. Existing laws and regulations are to be observed in all cases by our customers. This also applies in respect to any rights of third parties. Our information and advice do not relieve our customers of their own responsibility for checking the suitability of our products for the envisaged purpose.

The information in this document is subject to change without notice and should not be construed as a commitment by the manufacturing or selling entity, or an affiliate. We assume no responsibility for any errors that may appear in this document.

## Technical Assistance

Visit the tech service page at [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice).

## Terms and Conditions of Sale

Warranty, use restrictions, and other conditions of sale may be found at [SigmaAldrich.com/terms](https://SigmaAldrich.com/terms).

## Contact Information

For the location of the office nearest you, go to [SigmaAldrich.com/offices](https://SigmaAldrich.com/offices).

The life science business of Merck KGaA, Darmstadt, Germany  
operates as MilliporeSigma in the U.S. and Canada.

MilliporeSigma, and Sigma-Aldrich are trademarks of Merck KGaA, Darmstadt, Germany or its affiliates. All other trademarks are the property of their respective owners. Detailed information on trademarks is available via publicly accessible resources.

© 2019-2024 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved.