

03708 *Listeria* mono Enrichment Supplement I

For selective differentiation of *Listeria monocytogenes* from other *Listeria* species as per ISO Committee.

Composition

(per vial, sufficient for 500 ml medium)

L-phosphatidylinositol	1.0 g
Distilled water	25 ml

Directions

Thaw the contents of 1 vial of L. mono Enrichment Supplement I at room temperature. Aseptically add sterile contents of 1 vial of L. mono Enrichment Supplement I to 460 ml of sterile, molten L. mono Differential Agar Base (Cat. No. 77408) along with sterile rehydrated contents of L. mono Selective Supplement I (Cat. No. 92301) and L. mono Selective Supplement II (Cat. No. 91603). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

