


Product Information

COCAINE HYDROCHLORIDE

Methanol Solution


Product Number	C 1528	Synonyms	3-(Benzoyloxy)-8-azabicyclo-[3.2.1]octane-2-carboxylic acid methyl ester hydrochloride; l-Cocaine hydrochloride; Ecgonine methyl ester bezoate hydrochloride; Benzoylmethylecgonine hydrochloride
Lot Number	114K8804	Stability	If stored at 2-8°C in an air tight container, protected from light, decomposition of the product is less than 1% in 24 months. After opening, concentration may change due to loss of solvent.
CAS Number	[53-21-4]	Expiration Date	December 2007
Concentration	1.13 mg Cocaine Hydrochloride (equivalent to 1.01 mg free base) /mL GC grade methanol (conc. verified by UV)	Storage	2-8°C.
Molecular Formula	C ₁₇ H ₂₁ NO ₄ • HCl		
Molecular Weight	339.8		
DEA Class	Exempt preparation of Class II		

70 eV Electron Ionization Mass Spectrum



Hewlett-Packard 5890 Series II Gas Chromatograph
 Hewlett-Packard 5971 Mass Selective Detector
 Hewlett-Packard G1034C MS ChemStation (HP-DOS series)

FT-IR




Mattson Galaxy 3020 FT-IR
 KBr Pellet

COCAINE HYDROCHLORIDE

Methanol Solution


Product Number C 1528

Lot Number 114K8804



Ultraviolet Spectrum


Peak	233.8 nm
Absorbance	0.915
Solvent	0.2 N H ₂ SO ₄
Dilution	1:49 (v/v)



Bio-Tek Uvikon XL

Capillary Gas Chromatography

Column Supelco SPB-1 (Cat. No. 2-4162)
I.D. 0.20 mm **Length** 15 m **Film Thickness** 0.20 μm
Oven Temperature 50°C (2 min.), then 15°/min. to 300°C, hold 11 min.
Injector Temperature 250°C, Split
Detector Hewlett-Packard 5971
Mass Selective Detector (70 eV, 40-500 m/z)



Hewlett-Packard 5890

C 1528/114K8804/11/30/04

This product is not for *in vitro* diagnostic procedures. For R&D use only. Not for drug, household or other uses.

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications.

Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply.

Please see reverse side of the invoice or packing slip.