

## Product Information

**Anti-PSPC1 antibody, Mouse monoclonal**  
clone 1L4, purified from hybridoma cell culture

Catalog Number **SAB4200503**

### Product Description

Anti-PSPC1 (mouse IgG1 isotype) is derived from the hybridoma 1L4 produced by the fusion of mouse myeloma cells and splenocytes from BALB/c mice immunized with a synthetic peptide corresponding to a sequence at the C-terminal region of mouse PSPC1 (GeneID: 66645), conjugated to KLH. <sup>1</sup> The corresponding sequence is identical in rat and human PSPC1. The isotype is determined by ELISA using Mouse Monoclonal Antibody Isotyping Reagents, Catalog Number ISO2. The antibody is purified from culture supernatant of hybridoma cells grown in a bioreactor.

Anti-PSPC1 recognizes human and mouse PSPC1. The antibody can be used in several immunochemical techniques including immunoblotting (~60 kDa) and immunofluorescence.

PSPC1 (paraspeckle protein 1), also known as PSP1, is a member of the DBHS (Drosophila behavior, human splicing) protein family, consisting also of NONO (non-POU-domain-containing, octamer binding protein) and SFPQ (splicing factor, praline- and glutamine-rich). NONO, PSPC1 and SFPQ are the core protein components of paraspeckles. DBHS proteins bind both double- and single-stranded DNA and RNA, and are involved in various aspects of RNA and DNA metabolism, such as transcription, pre-mRNA 3' processing, transcription termination and mRNA splicing. These multi-functional nuclear proteins, under normal conditions, cycle between the nucleoplasm, paraspeckles and the nucleolus. DBHS proteins are involved in several biological processes, including the regulation of circadian rhythm, carcinogenesis and progression of cancer. PSPC1 activates androgen receptor mediated transcription. PSPC1 forms a heterodimer with NONO that localizes to paraspeckles in a RNA Polymerase II-dependent manner.<sup>1-5</sup>

### Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide.

©2021 Sigma-Aldrich Co. LLC. All rights reserved. SIGMA-ALDRICH is a trademark of Sigma-Aldrich Co. LLC, registered in the US and other countries. Sigma brand products are sold through Sigma-Aldrich, Inc. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see product information on the Sigma-Aldrich website at [www.sigmaaldrich.com](http://www.sigmaaldrich.com) and/or on the reverse side of the invoice or packing slip.

Antibody Concentration: ~ 1.0 mg/mL

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

### Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze at -20 °C in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

### Product Profile

Immunoblotting: a working concentration of 0.5-1 µg/mL is recommended using whole extracts of mouse Hepa1-6 cells.

Immunofluorescence: a working concentration of 5-10 µg/mL is recommended using human HeLa cells.

**Note:** In order to obtain the best results using various techniques and preparations, we recommend determining optimal working dilutions by titration.

### References

1. Kuwahara, S., et al., *Biol. Reprod.*, **75**, 352-359 (2006).
2. Bond, C.S. and Fox, A.H., *J. Cell Biol.*, **186**, 637-644 (2009)
3. Shav-Tal, Y., and Zipori, D., *FEBS Lett.*, **531**, 109-114 (2002).
4. Passon, D.M., et al., *Proc. Natl. Acad. Sci. USA*, **109**, 4846-4850 (2012).
5. Fox, A.H., et al., *Mol. Biol. Cell*, **16**, 5304-5315 (2005).

RC,ST,RC,PHC 04/21-1