

67706 Campylobacter Millichrome™ plus Supplement for 25L

Selective supplement for detection of thermotolerant Campylobacter.

Package information:

(1 vial is sufficient for 25 L of final media)

Chromogenic and selective mix 5.25g

Directions

- In a transparent vessel, add 210 mg of CHROMagar™ Campylobacter Supplement in 10 ml of purified water.
- Swirl well until complete dissolution.
- Filter to sterilize at 0.45 µm.
- Add the 10 ml of the supplement solution to 1 liter of the melted CHROMagar™ Campylobacter base at 45-50 °C.
- Swirl or stir gently to homogenize.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

