3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

Product Information

ADAM-10, mouse recombinant, expressed in Sf21 cells

Catalog Number **A9975** Storage Temperature –20 °C

Synonyms: Kuzbanian, mammalian disintegrin metalloprotease, or myelin-associated metalloproteinase

Product Description

Recombinant Mouse ADAM-10 is produced from the DNA sequence encoding the human CD33 signal peptide and the ectodomain of recombinant mouse ADAM-10 (amino acid residues 19-673). It is expressed with a C-terminal histidine tag in *Sf* 21 insect cells. Purified recombinant mouse ADAM-10, predicted molecular mass of ~52 kDa and 469 amino acids, exists as the mature and active form with the N-terminal sequence of T₂₁₅TLAERNTCQ. By immunoblotting, the apparent molecular mass is ~60 kDa.

ADAM-10 is a member of the ADAM (a disintegrin and metalloprotease-like domain) family. ^{1,2} It contains a pro domain with a cysteine switch and furin cleavage sequence, catalytic domain with the zinc-binding site and Met-turn (expected for reprolysins), disintegrin-like domain, cystein-rich domain, EGF-like domain, transmembrane domain, and cytoplasmic domain. ADAM-10 is highly conserved. Mouse, rat, cattle, and human ADAM-10 share 97% amino acid identity. There is 45% identity between mouse and *Drosophila* ADAM-10.

ADAM-10 is widely expressed in tissues and can be found on the cell surface and within the cell. Active ADAM-10 has an important role in neurogenesis (processing notch, notch ligand delta, and amyloid protein precursor at the alpha site). ADAM-10 also cleaves myelin basic protein and type IV collagen.

Recombinant Mouse ADAM-10 is lyophilized from a 0.2 μ m filtered solution in 25 mM MES, 0.2 M sodium chloride, 5 mM cacium chloride, 2.5 mM zinc chloride, 0.005% Brij[®] L23, pH 5.0.

Purity: >90% (SDS-PAGE, visualized by silver stain)

Endotoxin level: ≤1 EU/μg (LAL [Limulus amebocyte lysate] method) Recombinant Mouse ADAM-10 cleaves the coumarin peptide substrate, (7-methoxycoumarin-4-yl) acetyl-Pro-Leu-Ala-Gln-Ala-Val-(3-[2, 4-dinitrophenyl]-L-2,3-diaminopropionyl)-Arg-Ser-Ser-Ser-Arg-NH₂, and is measured by fluorescence (excitation at 320 nm and emission at 405 nm). The activity is measured with 10 μ M coumarin substrate and 100 ng of enzyme (ADAM-10) in a total of 100 μ l reaction mixture (50 mM HEPES, 5 μ M zinc chloride, 0.01% Brij L23, pH 7.5) at 37 °C. Under these conditions the specific activity is >5 pmoles/min/ μ g. In order to achieve the optimal activity, it is strongly recommended to use the described assay conditions. Use of other conditions may result in lower activity.

Recombinant Mouse ADAM-10 releases TNF α from a fusion protein containing proTNF α . It also cleaves the oxidized insulin B-chain. The major cleavage sites are Ala¹⁴–Leu¹⁵ and Tyr¹⁶–Leu¹⁷.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

Preparation Instructions

Reconstitute the contents of the vial with 25 mM Tris, pH 7.5, to give the desired stock solution (0.1 mg/ml recommended).

Storage/Stability

Store the product at -20 °C. ADAM-10 remains active for greater than six months at -20 °C to -70 °C.

Upon reconstitution, this enzyme can be aliquoted and stored at -20 °C for up to one month and -70 °C for three months without detectable loss of activity. Repeated freezing and thawing is not recommended. Do not store in a "frost-free" freezer.

References

- 1. Pan, D., and Rubin, G.M., Cell, **90**, 271-280 (1997).
- 2. Rooke, J., et al., Science, 273, 1227-1231 (1996).
- 3. Fahrenholz, F., *et al.*, Ann. N.Y. Acad. Sci., **920**, 215-222 (2000).
- 4. Lammich, S., et al., Proc. Natl. Acad. Sci. USA, **96**, 3922-3927 (1999).
- 5. Millichip, M.I., *et al.*, Biochem. Biophys. Res. Comm., **245**, 594-598 (1998).
- Rosendahl, M.S., et al., J. Biol. Chem., 272, 24588-24593 (1997).

Brij is a registered trademark of Uniqema Americas, LLC.

JR,RC,CS,KAA,MAM 07/12-1