

75208 Egg Yolk Tellurite Emulsion

Sterile stabilized tellurite emulsion of egg yolk recommended for the identification of *Staphylococcus* species.

Composition

(100 ml per vial)

Egg Yolk 30 ml, Sterile Saline 64 ml, Sterile 3.5% Potassium tellurite Solution 6 ml; Store at 4°C

Directions

The bottles are to be well shaken to suspend any sediment. Heat the refrigerated Egg Yolk Tellurite Emulsion in a water bath to 45 to 50°C and aseptically add 50 ml in 950 ml of sterile, molten Baird Parker Agar (11705). Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

