

POLYCYSTIN-2 CONTROL PEPTIDE

CATALOG NUMBER:	AG565
LOT NUMBER:	
QUANTITY:	100 µg
CONCENTRATION:	2 mg/mL (100 µg/50 µL)
DESCRIPTION:	Polycystin-2 control peptide.
APPLICATIONS:	Antibody blocking (AB9088) Optimal working concentration must be determined by the end user.
PRESENTATION:	Liquid in sterile distilled water.
STORAGE/HANDLING:	Store at -20°C in undiluted aliquots for up to 6 months after date of receipt. Avoid repeated freeze/thaw cycles.

For research use only; not for use as a diagnostic.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µl or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.

© 2007: Millipore Corporation. All rights reserved. No part of these works may be reproduced in any form without permissions in writing.