


**Enzyme Assay of ARYLAMINE ACETYLTRANSFERASE
(EC 2.3.1.5)****PRINCIPLE:**

Abbreviations used:

Acetyl-CoA = Acetyl Coenzyme A

Dye = p-Nitroaniline

AAT = Arylamine Acetyltransferase

CoASH = Coenzyme A

CONDITIONS:

T = 25°C, pH = 8.0, $A_{400\text{nm}}$, Light path = 1 cm

METHODS: Continuous Spectrophotometric Rate Determination

REAGENTS:

- A. 100 mM Potassium Phosphate Buffer, pH 8.0 at 25°C
(Prepare 100 ml in deionized water using Potassium Phosphate, Monobasic, Anhydrous, Sigma Prod. No. P-5379. Adjust to pH 8.0 at 25°C with 1 M KOH.)
- B. 100 mM β -Mercaptoethanol Solution (β -ME)
(Prepare 5 ml in deionized water using 2-Mercaptoethanol, Sigma Prod. No. M-6250.)
- C. 30 mM Ethylenediaminetetraacetic Acid Solution (EDTA)
(Prepare 5 ml in deionized water using Ethylenediaminetetraacetic Acid, Disodium Salt, Dihydrate, Sigma Stock No. ED2SS.)
- D. 1.5 mM p-Nitroaniline Solution (Dye)
(Prepare 5 ml in deionized water using p-Nitroaniline, Sigma Prod. No. N-2128.)
- E. 6.0 mM Acetyl Coenzyme A Solution (Acetyl CoA)
(Prepare 1 ml in deionized water using Acetyl Coenzyme A (C2:0), Sodium Salt, Sigma Prod. No. A-2056. **PREPARE FRESH.**)

**Enzyme Assay of ARYLAMINE ACETYLTRANSFERASE
(EC 2.3.1.5)**

REAGENTS: (continued)

- F. Arylamine Acetyltransferase Enzyme Solution
(Immediately before use, prepare a solution containing 100 mg/ml of Arylamine Acetyltransferase in cold Reagent A.)

PROCEDURE:

Pipette (in milliliters) the following reagents into suitable containers:

	<u>Test</u>	<u>Blank</u>
Reagent A (Buffer)	1.50	1.50
Reagent B (β-ME)	0.15	0.15
Reagent C (EDTA)	0.10	0.10
Reagent D (Dye)	0.20	0.20
Reagent F (Enzyme Solution)	1.00	1.00

Mix by inversion and equilibrate to 25°C. Monitor the A_{400nm} until constant (approx. 4 minutes) using a suitably thermostatted spectrophotometer. Then add:

Reagent E (Acetyl CoA)	0.05	-----
Deionized Water	-----	0.05

Immediately mix by inversion and record the increase in A_{400nm} for approximately 5 minutes. Obtain the $\Delta A_{400nm}/\text{minute}$ using the maximum linear rate for both the Test and Blank.

CALCULATIONS:

$$\text{Units/ml enzyme} = \frac{(\Delta A_{400nm} \text{ Test} - \Delta A_{400nm} \text{ Blank})(1000)(3)(df)}{(11.19)(1)}$$

1000 = Conversion from μmoles to nanomoles

3 = Volume (in milliliters) of assay

df = Dilution factor

11.19 = Millimolar extinction coefficient of acetylated dye

1 = Volume (in milliliters) of enzyme used

$$\text{Units/mg solid} = \frac{\text{units/ml enzyme}}{\text{mg solid/ml enzyme}}$$

**Enzyme Assay of ARYLAMINE ACETYLTRANSFERASE
(EC 2.3.15)**

CALCULATIONS: (continued)

$$\text{Units/mg protein} = \frac{\text{units/ml enzyme}}{\text{mg protein/ml enzyme}}$$

UNIT DEFINITION:

One unit will acetylate 1.0 nanomole of p-nitroaniline per minute at pH 8 at 25°C.

FINAL CONCENTRATION:

In a 3.00 ml reaction mix, the final concentrations are 83 mM potassium phosphate, 5 mM β -mercaptoethanol, 1 mM ethylenediaminetetraacetic acid, 0.1 mM p-nitroaniline, 0.1 mM acetyl coenzyme A and 100 mg of arylamine acetyltransferase.

REAGENTS:

Brooks, S.P.J. and Storey, K.B., (1993) *Analytical Biochemistry* **212**, 452-456

NOTES:

1. This assay is based on the cited reference.
2. Where Sigma Product or Stock numbers are specified, equivalent reagents may be substituted.

Sigma warrants that the above procedure information is currently utilized at Sigma and that all Sigma-Aldrich, Inc. products conform to the information in this and other Sigma-Aldrich, Inc. publications. Purchaser must determine the suitability of the information and product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.