



## QUALITY CONTROL RANGES

### MILLIPLEX® MAP Human Neurodegenerative Disease Magnetic Bead Panel 1 Catalog # HNDG1MAG-36K

#### Control Catalog # HNDG1-6036-1 Control Lot # HNDG1-110.210

| <u>Hormone</u> | <u>QC Level</u> | <u>Lot #</u> | <u>Expected Range</u> | <u>Units</u> |
|----------------|-----------------|--------------|-----------------------|--------------|
| APOA-1         | Control I       | HNDG1-110    | 12 - 25               | ng/mL        |
|                | Control II      | HNDG1-210    | 61 - 126              | ng/mL        |
| APOC-3         | Control I       | HNDG1-110    | 0.19 - 0.39           | ng/mL        |
|                | Control II      | HNDG1-210    | 0.9 - 1.8             | ng/mL        |
| APOE           | Control I       | HNDG1-110    | 0.61 - 1.7            | ng/mL        |
|                | Control II      | HNDG1-210    | 3.8 - 8.0             | ng/mL        |
| HCFH           | Control I       | HNDG1-110    | 5.4 - 11              | ng/mL        |
|                | Control II      | HNDG1-210    | 25 - 52               | ng/mL        |
| HCC3           | Control I       | HNDG1-110    | 0.80 - 1.7            | ng/mL        |
|                | Control II      | HNDG1-210    | 3.9 - 8.2             | ng/mL        |
| HA2MG          | Control I       | HNDG1-110    | 7.7 - 16              | ng/mL        |
|                | Control II      | HNDG1-210    | 37 - 76               | ng/mL        |

HNDG1-6036-1: HNDG1-110.210 (26 June 2024)

EMD Millipore • 1-800-645-5476 • 1-800-645-5439 (FAX) • [www.emdmillipore.com](http://www.emdmillipore.com)

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA