



QUALITY CONTROL RANGES

HUMAN FGF-21 ELISA CATALOG # EZHFGF21-19K

HUMAN FGF-21 STANDARD, LYOPHILIZED CATALOG# E8019-K

STANDARD	HUMAN FGF-21 VALUE
Lot Number. 8019-12K	2000 pg/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6019-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	FGF21-112	82 – 170	pg/mL
QC 2	FGF21-212	532 – 1105	pg/mL

Cat# E6019-K: FGF21-112.212 (rev. 11 Oct. 2023)
EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX)

• www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.
© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA