



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone 800-325-5832 • (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

Monoclonal Anti-hPlk1 (Polo Like Kinase 1)

Clone 36-298

produced in mouse, purified immunoglobulin

Catalog Number **P6123**

Product Description

Monoclonal Anti-hPlk1 (Polo Like Kinase 1) (mouse IgG1 isotype) is derived from the 36-298 hybridoma produced by the fusion of mouse myeloma cells (PA1) and splenocytes from BALB/c mice immunized with human full length Plk1. The isotype is determined using Sigma ImmunoType™ Kit, Catalog Number ISO1, and by a double diffusion immunoassay using Mouse Monoclonal Antibody Isotyping Reagents, Catalog Number ISO2.

Monoclonal Anti-hPlk1 recognizes human, mouse and rat Plk1 (approx. 63 kDa). The product is useful in ELISA, immunoblotting, immunoprecipitation and immunocytochemistry (paraformaldehyde fixation). This antibody, clone 36-298, is recommended over clone 35-206, Catalog Number P5998, for immunocytochemistry assays.

Polo-like kinases (Plks) are a family of serine/threonine kinases that regulate cellular events in mitosis and cytokinesis. The Plk family is structurally and functionally related to the polo gene in *Drosophila*.^{1,2} Four Plks were found to be expressed in mammals: Plk1, Plk2 or Snk, Plk3 or Prk/Fnk and Sak.^{2,3} Plk1 activity is regulated by its phosphorylation and normally appears during mitosis. Plk1 localizes during mitosis in centrosomes, kinetochores and in the central spindle. Plk1 phosphorylates a number of substrate proteins such as Cdc25C and cyclin B1. In addition to its role in regulation of cell cycle, it was suggested that Plk1 has a role in tumor progression. Studies have shown that Plk1 is a target of the DNA damage checkpoint on the basis that its activity is inhibited by DNA damage in G₂ and mitosis. Moreover, Plk1 was found to bind to the sequence-specific DNA-binding domain of p53. This interaction between p53 and Plk1 inhibits the transcriptional activation and the pro-apoptotic functions of p53.^{4,6}

Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

Antibody concentration: ~2.0 mg/mL

Precautions and Disclaimer

Due to the sodium azide content a material safety sheet (MSDS) for this product has been sent to the attention of the safety officer of your institution.

Consult the MSDS for information regarding hazardous and safe handling practices.

Storage/Stability

For continuous use, store at 2-8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing, or storage in "frost-free" freezers, is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

Product Profile

A working concentration of 0.5-1.0 µg/mL is determined by immunoblotting, using a whole extract of HeLa cell line (Human cervix epitheloid carcinoma).

Note: In order to obtain best results in different techniques and preparations we recommend determining optimal working concentration by titration test.

References

1. Glover, D.M., et al., *J. Cell Biol.*, **135**, 1681-1684 (1996).
2. Cheng, K-Y., et al., *EMBO J.*, **22**, 5757-5768 (2003).

3. Van Vugt, M.A.T.M, et al., *J. Biol. Chem.*, **279**, 36841-36854 (2004).
4. Neef, R., et al., *J. Cell Biol.*, **162**, 863-875 (2003).
5. Smits, V.A.J., et al., *Nat. Cell Biol.*, **2**, 672-676 (2000).
6. Ando, K., et al., *J. Biol. Chem.*, **279**, 25549-25561 (2004).

RF,AH,PHC 03/06-1

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.