



QUALITY CONTROL RANGES

SARS-COV2 ANTI-RBD IgG ELISA KIT CATALOG # EZRBDG-110K

STANDARD, LYOPHILIZED CATALOG # E8110-K

STANDARD	SARS-COV2 ANTI-RBD IgG VALUE
LOT NUMBER: 8110-1	100 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **1.0mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG #E6110-K

Reconstitute each vial with **1.0mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	6110-101	5.3 – 11	ng/mL
QC 2	6110-201	18 - 37	ng/mL

E6110-K: 6110-101.201 (31 March 2021)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com