



MOUSE ANTI-INFLUENZA B MONOCLONAL ANTIBODY

CATALOG NUMBER:	MAB8661
LOT NUMBER:	2216879
QUANTITY:	100 µL
SPECIFICITY:	All subtypes of Influenza B. Blend of two clones that recognize the Influenza B nucleoprotein.
IMMUNOGEN:	Influenza Blend.
ISOTYPE:	IgG ₁
APPLICATIONS:	Indirect immunofluorescence at 1:100. Optimal working dilutions must be determined by end user.
SPECIES REACTIVITIES:	No cross-reactivity to influenza A, parainfluenza 1, 2 or 3, adenovirus or RSV.
FORMAT:	Ascites fluid.
PRESENTATION:	Liquid containing 0.1% sodium azide.
STORAGE/HANDLING:	Maintain at -20°C in convenient undiluted aliquots for up to 12 months after date of receipt. Avoid repeated freeze/thaw cycles.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

**FOR RESEARCH USE ONLY; NOT FOR USE IN DIAGNOSTIC
PROCEDURES. NOT FOR HUMAN OR ANIMAL CONSUMPTION**

Unless otherwise stated in our catalog or other company documentation accompanying the product(s), our products are intended for research use only and are not to be used for any other purpose, which includes but is not limited to, unauthorized commercial uses, in vitro diagnostic uses, ex vivo or in vivo therapeutic uses or any type of consumption or application to humans or animals.

©2002 2012 EMD Millipore . All rights reserved. No part of these works may be reproduced in any form without permission in writing.

EMD Millipore Corporation, 28820 Single Oak Drive • Temecula, CA 92590 USA 1-800-437-7500
Technical Support: T: 1-800-MILLIPORE (1-800-645-5476) • F: 1-800-437-7502
www.millipore.com

December 2012
MAB8661CA; Rev. A