



QUALITY CONTROL RANGES

Rat/Mouse Growth Hormone ELISA Catalog #: EZRMGH-45K

STANDARD, LYOPHILIZED CATALOG #: E8045-K

STANDARD	RMNPY VALUE
LOT NUMBER: 845-18K	57 ng/mL

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG # E6045-K

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	RMGH-118	1.6 – 3.4	ng/mL
QC 2	RMGH-218	12 – 25	ng/mL

E6045-K: RMGH-118.218 (Rev 10 May 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA