

**MOUSE ANTI-RESPIRATORY SYNCYTIAL VIRUS (RSV)
MONOCLONAL ANTIBODY**

CATALOG NUMBER: MAB858-2

LOT NUMBER: xxxx

QUANTITY: 100 µg

CONCENTRATION: 1 mg/mL

SPECIFICITY: Glycoprotein, 89 kD. Recognizes both A and B strains.

IMMUNOGEN: A2

ISOTYPE: IgG_{1,k}

CLONE NAME: 131/2G

APPLICATIONS: ELISA
Indirect Immunofluorescence (fresh frozen tissue sections)
Final working dilutions must be determined by end user.

FORMAT: Purified immunoglobulin at 1 mg/mL in 0.02M PB with 0.25M NaCl, pH=7.6 containing 0.1% sodium azide.

PRESENTATION: Liquid.

STORAGE/HANDLING: Maintain at 2-8° C in aliquots for up to 12 months.
Do not store in diluted format.

Important Note: *During shipment, small volumes of product will occasionally become entrapped in the seal of the product vial. For products with volumes of 200 µL or less, we recommend gently tapping the vial on a hard surface or briefly centrifuging the vial in a tabletop centrifuge to dislodge any liquid in the container's cap.*

For research use only; not for use as a diagnostic.