



Quality Control Ranges
Milliplex PLEXpedition Magnetic Bead Screening Panel
Catalog No. HPLX1-115SP, HPLX1-115SP-PX, HPLX1-115SP-PX80
Control Catalog No. HPLX1-6010-06 Lot No. PLX16-101 and PLX16-201

<u>Bead</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>	<u>Bead</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HBNP	Control 1	495 - 1027	pg/mL	HOPG	Control 1	129 - 268	pg/mL
	Control 2	2517 - 5228	pg/mL		Control 2	638 - 1326	pg/mL
HCKMB	Control 1	3975 - 8256	pg/mL	HTRPNI	Control1	85 - 177	pg/ml
	Control 2	20637 - 42861	pg/mL		Control 2	456 - 947	pg/ml
HDKK1	Control 1	46 - 96	pg/mL	HNTPBNP	Control1	62 - 130	pg/ml
	Control 2	237 - 492	pg/mL		Control 2	314 - 652	pg/ml
HACTH	Control 1	31 - 64	pg/mL	HOC	Control1	493 - 1025	pg/ml
	Control 2	152 - 317	pg/mL		Control 2	2549 - 5293	pg/ml
HFABP3	Control 1	299 - 620	pg/mL	HOSTNCTN	Control1	4764 - 9894	pg/ml
	Control 2	1579 - 3279	pg/mL		Control 2	22937 - 47638	pg/ml
HFABP4	Control 1	81 - 169	pg/mL	HSOST	Control1	615 - 1277	pg/ml
	Control 2	435 - 904	pg/mL		Control 2	3153 - 6549	pg/ml
HOPN	Control 1	493 - 1024	pg/mL				
	Control 2	2540 - 5276	pg/mL				

HPLX1-6010-06: PLX16-101.201 (04 September 2025) • [SigmaAldrich.com/support](https://sigmaaldrich.com/support) • [SigmaAldrich.com](https://sigmaaldrich.com)