


30 Liter Polyethylene Tank

(Cat. No. TANKPE030)

General Description

Purifying water is useless if the quality of the pure water degrades during storage. Merck's 30 L polyethylene (PE) tanks are designed to maintain consistent purity of the stored water and provide effective protection against airborne contaminants. A Millipak® vent filter can be installed on the tank to avoid contamination due to the air entering the tank.



Installation requirements

30 L tank mounting dimensions: 697 mm (27.5 in.) H x 382 mm (15 in.)

Installation options

Wall-mounting options for a single unit require a wall capable of supporting 30 kg (66. 14 lb.).

A piece of plywood 2 ft. wide x 2 ft. $\frac{3}{4}$ in. high is recommended.

Wall Mounting Bracket (TANKFIX01) is required.


Water connection 5 $\frac{1}{4}$ in. FNPT (threaded plugs not included).


Electrical specifications

None

Operating weight

Maximum usable capacity	25 L
Weight full	30 kg (66.14 lb.)





For More Information

Please visit our website

[MerckMillipore.com/labwater](https://www.MerckMillipore.com/labwater)

Merck KGaA
Frankfurter Strasse 250
64293 Darmstadt
Germany

