



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information

LAMBDA PHAGE DNA

Product No. **D 3779**

Lot No. 015k6009

Storage: Below 0 °C

PRODUCT SUMMARY

Source: *Escherichia coli* W3110

Buffer: 10 mM Tris-HCl, pH 8.0, 1 mM EDTA

Concentration: 506 µg/ml or 10.1 A_{260} units/ml

Ratio: A_{260}/A_{280} : 1.8

M.W.: 31.5×10^6 Daltons

BASE PAIRS: 48,502

Unique Sites: Apa I, Nae I, Nar I, Nhe I, PaeR7 I, SnaB I,
Xba I, Xho I

COMMENTS

Strain $\lambda c1857$, Sam 7

Note: One A_{260} unit is equivalent to approximately 50 µg of DNA.

Reference:

Sanger, F., et al., J. Mol. Biol., **162**, 729 (1982).

3/05