



3050 Spruce Street
Saint Louis, Missouri 63103 USA
Telephone (800) 325-5832 (314) 771-5765
Fax (314) 286-7828
email: techserv@sial.com
sigma-aldrich.com

Product Information


DBO-83

Product Number **D 7938**

Storage Temperature: Room Temperature

CAS#: 195211-53-1

Synonym: 3,8-Diazabicyclo[3.2.1]octane,3-(6-chloro-3-pyridazinyl)-,dichloride



2HCl

Product Description

Molecular Formula: $C_{10}H_{13}ClN_4 \cdot 2HCl$

Molecular Weight: 297.69 (anhyd.)

DBO-83 is a novel nicotinic acetylcholine receptor agonist, structurally related to epibatidine, with antinociceptive properties in rodents.

Preparation Instructions

Soluble in water (12 mg/ml).

Storage/Stability

Store tightly sealed at room temperature.

References

1. Lloyd, G.K. et al., J. Pharmacol. Exp. Ther., **292**, 461-67 (2000).
2. Ghelardini, C. et al., Drug Dev. Res., **40**, 251-258 (1997).

mje 8/00

Sigma brand products are sold through Sigma-Aldrich, Inc.

Sigma-Aldrich, Inc. warrants that its products conform to the information contained in this and other Sigma-Aldrich publications. Purchaser must determine the suitability of the product(s) for their particular use. Additional terms and conditions may apply. Please see reverse side of the invoice or packing slip.