3050 Spruce Street, St. Louis, MO 63103 USA
Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757
email: techservice@sial.com sigma-aldrich.com

## **Product Information**

# Puromycin dihydrochloride from Streptomyces alboniger

Catalog Number **P7255** Storage Temperature –20 °C

CAS RN: 58-58-2 (dihydrochloride) 53-79-2 (free base)

Synonyms: Stylomycin dihydrochloride;

(S)-3'-[[2-amino-3-(4-methoxyphenyl)-1-oxopropyl]-amino]-3'-deoxy-N,N-dimethyl-adenosine

dihydrochloride

Molecular formula:  $C_{22}H_{29}N_7O_5 \cdot 2HCI$ 

Molecular weight: 544.43

Melting Point:  $^{1,2}$  175.5–177 °C (free base) pK<sub>a</sub> values:  $^{2,3}$  6.8, 7.2 Extinction coefficient:  $^2$  20.3 (275 nm, 0.1 N NaOH) Extinction coefficient: 19.5 (268 nm, 0.1 N HCI)  $\left[\alpha\right]^{25}$ :  $^4$  –11° (ethanol)

#### **Product Description**

Puromycin dihydrochloride belongs to the aminonucleoside family of antibiotics and is isolated from *Streptomyces alboniger*. Since the partial structure of this antibiotic showed it to be a purine derivative, puromycin was assigned as its generic name.

Puromycin is a broad spectrum antibiotic and antibacterial agent. It is active against Gram-positive microorganisms, less active against acid-fast bacilli, and weakly active against Gram-negative microorganisms. <sup>2.5</sup> It prevents growth of bacteria, protozoa, algae, and mammalian cells. Puromycin also possesses antitrypanosomal activity. It acts very quickly and can kill 99% of the cells within 2 days. It also exhibits antitumor activity in studies on brain tumor cells. <sup>6</sup>

Puromycin is a protein synthesis inhibitor that causes premature chain termination by acting as an analog of the 3'-terminal end of aminoacyl-tRNA. It has been found to inhibit the incorporation of <sup>14</sup>C-L-leucine into protein in a cell-free preparation from rat liver. It was postulated the inhibition was due to the fact that no <sup>14</sup>C-L-leucine was transferred from soluble (or transfer) ribonucleic acid—<sup>14</sup>C-L-leucine to microsomal protein. It has been used to study transcriptional regulatory mechanisms that control the sequential and coordinate expression of genes during cell differentiation. <sup>7,8</sup>

The optimal working concentration of puromycin varies between cell lines. Usually, the working concentration for eukaryote cell culture is 1–10 µg/ml.

#### **Precautions and Disclaimer**

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

#### **Preparation Instructions**

The product is soluble in water (50 mg/ml), yielding a clear, colorless to faint yellow solution. The stock solution may be passed through a 0.22  $\mu$ m filter and stored in aliquots at –20 °C.

It is also soluble in methanol (10 mg/ml).

### Storage/Stability

Store the product at  $-20~^{\circ}$ C. Under these conditions, the powder remains active for four years.

#### References

- 1. Fryth, P.W., et al., *J. Amer. Chem. Soc.*, **80**, 2736-2740 (1958).
- 2. Data for Biochemical Research, 3rd ed., 280-281 (1986).
- 3. Waller, C.W., et al., *J. Amer. Chem. Soc.*, **75**, 2025 (1953).
- 4. Korzybski, T., et al., *Antibiotics*, **II**, 1173-1180 (1978).
- Antibiotics: origin, nature, and properties, Korzybski, T., et al., American Society for Microbiology (Washington, DC: 1978), pp. 1173-1180.
- 6. Lee, Y.S., and Wurster, R.D., *Cancer Lett.*, **93**, 157-163 (1995).
- 7. Manabe I. And Owens, G.K., *Circ. Res.*, **88**, 1127-34 (2001).
- 8. Marchetti S., et al., *J. Cell Sci.*, **115**, 2075-85 (2002).

ES-S,PHC,MAM 01/09-1