

**QUALITY CONTROL RANGES  
CANINE C-PEPTIDE RIA  
CATALOG # CCP-24HK**



**CALIBRATOR, LYOPHILIZED  
CATALOG #8024-K**

CALIBRATOR	CANINE C-PEPTIDE VALUE
LOT NUMBER: CCP80-13K	22 ng/mL

The exact value of lyophilized calibrator is shown above.  
Reconstitute with **2 mL** distilled or deionized water.

**QUALITY CONTROLS 1 AND 2, LYOPHILIZED  
CATALOG #6000L-K**

Reconstitute each vial with **1 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGES	UNITS
QC 1	UV-128	0.91 – 1.9	ng/mL
QC 2	UV-228	3.4 – 7.1	ng/mL

6000L-K: UV-128.228 (29 January 2019)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) [www.emdmillipore.com](http://www.emdmillipore.com)