



QUALITY CONTROL RANGES

**MILLIPLEX[®] MAP Human TIMP Panel 1 Magnetic
Kit Catalog # HTMP1MAG-54K**

**MILLIPLEX[®] MAP Human TIMP Panel 2 Magnetic
Kit Catalog # HTMP2MAG-54K**

**Control Catalog # HTIMP-6054
Lot# HTMP-114 and HTMP-214**

<u>Hormone</u>	<u>QC Level</u>	<u>Expected Range</u>	<u>Units</u>
HTIMP-1	Control I	257 - 534	pg/mL
	Control II	1288 - 2675	pg/mL
HTIMP-2	Control I	634 - 1317	pg/mL
	Control II	3206 - 6658	pg/mL
HTIMP-3	Control I	1213 - 2520	pg/mL
	Control II	6582 - 13670	pg/mL
HTIMP-4	Control I	138 - 287	pg/mL
	Control II	727 - 1509	pg/mL

Note: The Quality Control Ranges are generated with overnight assay format using assay buffer provided in the kit. Quality Control values in culture media are not tested.

HTIMP-6054: HTMP-114.214 (15 August 2023)

EMD MILLIPORE • 1-800-645-5476 • 1-800-645-5439 (FAX) • www.emdmillipore.com

The M logo is a trademark of Merck KGaA, Darmstadt, Germany.

© 2013 EMD Millipore Corporation, Billerica, MA 01821 USA