


# phosphoramidites and reagents

## for Oligonucleotide Synthesis

Serving Your Needs across the Entire  
Oligonucleotide Synthesis Workflow



# Proligo® Reagents for Oligonucleotide Synthesis

Let us support your entire oligonucleotide production workflow, from synthesis, cleavage and deprotection through to desalting, purification and QC analysis. Many of our oligo synthesis columns and reagents are specifically configured to work with your synthesizer (ABI 3900, 39X, Expedite, MerMade, AKTA® and OligoPilot®) from R&D to production.

## The Oligonucleotide Synthesis Cycle

### **We offer a comprehensive range of high-quality products and services for oligonucleotide synthesis:**


- Our products span the synthesis workflow, from standard and specialty phosphoramidites, to solid supports, columns, linkers and liquid formulations.
- Products are available to fit directly on the instruments with specialized bottles, or in bulk packaging to feed larger scale processes.
- We offer custom manufacturing, packaging, and global distribution capabilities.

### **Use our technical and regulatory know-how, global quality standards, and supply chain expertise to help you:**


- Improve efficiency
- Bring your products to market faster
- Optimize value

To view our products used in each stage of the synthesis cycle, visit [SigmaAldrich.com/proligo](https://SigmaAldrich.com/proligo)


### Cleavage and Deprotection after x cycle




## CPG loaded with nucleoside





Cap A and B solution




### Deblock solution




Oxidizer solution



Cap A and B solution



## DNA/RNA Phosphoramidites Labels/Dyes/Linkers



\* = B<sub>r</sub> is equivalent to BL in Cycle

## Synthesis

Our complete product portfolio for solid-phase oligonucleotide synthesis is comprised of high-quality standard and specialty phosphoramidites, solid supports, columns, linkers and liquid formulations.

### Raw Materials for Oligonucleotide Synthesis

- Controlled Pore Glass (CPG) in bulk or packed Columns for all major synthesizers
- DNA and RNA Phosphoramidites in standard and fast deprotection as well as custom solutions
- Linkers, Dyes & Labels
- Liquid reagents, from bottles to bulk, for detritylation, activation, coupling, oxidation and capping

### Cleavage and Deprotection

Use our high-quality solvents and reagents for cleavage and deprotection to help you improve your yields and reduce unwanted side reactions.

- Ethylenediamine
- Methylamine
- Ammonium Hydroxide
- Tetrabutylammonium Bromide
- Hydrogen Fluoride Pyridine
- Triethylamine

### Desalting and Purification

We can help with the removal of organic salts through a variety of methods, such as:

- Chromatography columns
- Nuclease-free molecular biology reagents
- Gel electrophoresis
- Organic solvents for extractions
- Ethanol for precipitation

## Desalting and Purification

Oligonucleotides are best purified by different methods or a combination of techniques, depending on the intended use.

## Analysis

Materials to analyze composition and purity.

### HPLC

- Reversed Phase Columns
- Ion Exchange Columns
- High Purity Solvents

### Analytical Consumables for:

- HPLC, UHPLC, LC-MS
- GC, headspace GC, GC-MS
- TLC, HPTLC
- Karl Fischer titration
- NMR
- UV-Vis
- FTIR

For more on analytical applications, visit [SigmaAldrich.com/brighter](https://SigmaAldrich.com/brighter)

## A Complete Line of Solvents and Labware for Your Oligonucleotide Synthesis

### Solvents

We have an expansive range of products, grades, and packaging options. Our state-of-the-art manufacturing capabilities provide reliable quality and lot-to-lot consistency. Products include:

- |                     |                     |
|---------------------|---------------------|
| • Acetone           | • Acetonitrile      |
| • Dichloromethane   | • Chloroform        |
| • Dimethylformamide | • Dimethylsulfoxide |
| • Diethyl Ether     | • Ethanol           |
| • Ethyl Acetate     | • Hexane            |
| • Methanol          | • 2-Propanol        |
| • Tetrahydrofuran   | • Toluene           |

For more on solvents, visit [SigmaAldrich.com/solvents](https://SigmaAldrich.com/solvents)

### Labware

Over 22,000 products from the most-recognized brands in the world. We carry all the products that your lab needs. Product categories specific to the oligonucleotide market include:

- |                         |                               |
|-------------------------|-------------------------------|
| • PCR Supplies          | • Microtiter/Microwell Plates |
| • Centrifugation        | • Filtration                  |
| • Cell Culture Supplies | • Electrophoresis             |


For more on labware, visit [SigmaAldrich.com/labware](https://SigmaAldrich.com/labware)



# Your One-Stop-Shop for Lab Consumables

We are more than just your preferred supplier for products and services around the Oligonucleotide synthesis cycle, we can also support your entire workflow with our complete portfolio:

- Analytical products and tools for use in determining product composition and impurity profile
- Reagents and solvents for synthesis, purification and analytical purposes
- Labware from the most-recognized brands in the world



## Lab Consumables Portfolio

| Description                                     |   |
|---|---|
| <b>Analytical Reagents and Equipment</b>        |   |
| Reagents  | Test Kits and Indicators                            |
|   | Spectroscopy  |
|   | Derivatization                                      |
|   | Electrochem, Electrodes                             |
|   | Chromatography                                      |
| Electrochemistry, Electrodes, Sensors           |   |
| Titration and Karl Fischer Analysis             | Buffers, Indicators, Volumetric Solutions           |
| Lab Equipment                                   | Syringes, Vials                                     |
| Standards                                       | Certified Reference Materials                       |
|   | Physical Properties                                 |
|   | Chromatography and Spectroscopy                     |
| General Lab Equipment (Containers, Bottles)     |   |
| Genomics  | Nucleosides, Substrates                             |
| Cell Culture                                    | Reagents and Components                             |
| Biomolecules                                    |   |
| Proteins and Enzymes                            |   |
| <b>Chemical Synthesis</b>                       |   |
| Building Blocks                                 |   |
| Catalysts and Inorganics                        |   |
| Synthetic Reagents                              |   |
| <b>Lab Equipment</b>                            |   |
| Corning Plastic and Glassware, Centrifuge tubes |   |
| General Laboratory, Filtration                  |   |
| <b>PCR Supplies</b>                             |   |
| Centrifugation                                  |   |
| Cell Culture Supplies                           |   |
| Microtiter/Microwell Plates                     |   |
| Filtration                                      |   |
| <b>Lab Essentials</b>                           |   |
| Labbasics                                       | Inorganics Reagents                                 |
|   | Organics  |
|   | Acid/Bases  |
|   | Solvents  |
| Biobasics                                       | Buffers, Carbohydrates                              |
|   | Molecular Biology Reagents, Nucleosides/Nucleotides |
|   | Detergents  |

## To place an order or receive technical assistance

Phone: +49 40 79 702 250  
E-mail: [DEUHAM-CS@sial.com](mailto:DEUHAM-CS@sial.com)

Order/Customer Service: [SigmaAldrich.com/order](https://SigmaAldrich.com/order)  
Technical Service: [SigmaAldrich.com/techservice](https://SigmaAldrich.com/techservice)  
Safety-related Information: [SigmaAldrich.com/safetycenter](https://SigmaAldrich.com/safetycenter)

[SigmaAldrich.com/proligo](https://SigmaAldrich.com/proligo)

