

Product Information

Anti-Human Serum antibody produced in goat whole antiserum

Catalog Number **H9640**

Product Description

Anti-Human Serum was developed in goat using pooled normal human serum as the immunogen. The antiserum has been treated to remove lipoproteins.

Strong reactivity with normal human serum has been determined by immunoelectrophoresis (IEP). This antiserum has not been assayed for interspecies cross reactivity.

Reagents

Goat Anti-Human Serum is provided as a liquid containing 0.1% sodium azide as preservative.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

This goat antiserum was maintained at pH 5.0 for 40 minutes to meet U.S.D.A. requirements.

Storage/Stability

For continuous use, store at 2–8 °C for up to one month. For extended storage, the solution may be frozen in working aliquots. Repeated freezing and thawing is **not** recommended. Storage in "frost-free" freezers is **not** recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use.

DS,PCS,KMR,MAM 12/18-1