

Product No. I-4895
Lot 044H0298

Anti-Mouse Interleukin-9
Developed in Goat
IgG Fraction of Antiserum

Anti-Mouse Interleukin-9 is developed in goat using recombinant, mouse interleukin-9 (rmIL-9), expressed in the insect cell line *Sf* 21, as the immunogen. The product is purified by Protein G affinity chromatography. Goat Anti-Mouse IL-9 is provided lyophilized from phosphate buffered saline, pH 7.4, to which no preservatives have been added.

Description

Interleukin-9 is a T cell-derived¹ glycoprotein having multiple functions on cells of the lymphoid, myeloid, and mast cell lineages.² Anti-Mouse IL-9 neutralizes the bioactivity of rmIL-9 but not that of human, recombinant IL-9. In indirect ELISA and immunoblotting, this antibody displays approximately 25% cross-reactivity with human, recombinant IL-9. The product shows no cross-reactivity with other tested cytokines using indirect ELISA.

Performance

Anti-Mouse IL-9 was tested for its ability to neutralize the bioactivity of rmIL-9 in a cell proliferation assay using the factor-dependent cell line TS1.³ The ND₅₀ of the antibody is defined as the concentration of antibody resulting in a one-half maximal inhibition of bioactivity of rmIL-9, which is present at five times its own EC₅₀ (the concentration of rmIL-9 producing a one-half maximal bioactivity without antibody). In this bioassay, rmIL-9 was pre-incubated with various dilutions of the antibody for 1 hour at 37°C in a 96-well microtiter plate. Then, TS1 cells were added to each well to give a final concentration of 2.5 x 10³ cells/ml in 0.2 ml containing 0.05 ng/ml rmIL-9. This was incubated for 72 hours at 37°C in a 5% CO₂ humidified incubator and then pulsed for 4 hours with ³H-thymidine. Cells were harvested onto glass filters and the ³H-thymidine incorporation into DNA was measured.

Product Information

Mass/vial:	1 mg
Immunogen:	Mouse, recombinant IL-9 (rmIL-9)
Host animal:	Goat
Formulation:	Lyophilized from PBS without additives
Endotoxin:	<10 ng/vial by LAL method
Bioactivity:	ND ₅₀ = 4.1 µg/ml
Indirect ELISA:	0.5 µg/ml detects 0.15 ng/well of rmIL-9
Indirect	
Immunoblotting:	1 µg/ml antibody detects rmIL-9 at 1 ng/lane under both reducing and non-reducing conditions
Sterility:	0.2 µm-filtered, aseptic fill

Reconstitution and Use

To one vial of lyophilized powder, add 1 ml of sterile-filtered PBS to produce a 1 mg/ml stock solution of Anti-Mouse IL-9. If aseptic technique is used, no further filtration should be needed for use in cell culture environments.

Storage

Prior to reconstitution, store at -20°C for 6 months. Reconstituted product may be stored at 0-5°C for up to one month. For prolonged storage, freeze in working aliquots at -20°C. Avoid repeated freezing and thawing.

References

1. Yang, Y., et al., Blood, **74**, 1880 (1989).
2. Renauld, J., et al., Int. Rev. Exp. Path., **34A**, 99 (1993).
3. Uyttenhove, C., et al., Proc. Natl. Acad. Sci. USA, **85**, 6934 (1988).

Prepared for Sigma.

Sigma warrants that its products conform to the information contained in this and other Sigma publications. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.

Issued 06/94.