



QUALITY CONTROL RANGES

GLP-1 TOTAL ELISA CATALOG# EZGLP1T-36K, EZGLP1T-36BK

STANDARD, LYOPHILIZED CATALOG# E8036-K

STANDARD	GLP-1 TOTAL VALUE
Lot Number. 836-18K	1023 pM

The exact value of lyophilized standard is shown above.

Reconstitute with **0.5 mL** distilled or deionized water.

QUALITY CONTROLS 1 AND 2, LYOPHILIZED CATALOG# E6036-K

Reconstitute each vial with **0.5 mL** distilled or deionized water.

CONTROL	LOT NUMBER	EXPECTED RANGE	UNITS
QC 1	GLP1-118	25 - 52	pM
QC 2	GLP1-218	140 - 291	pM

E6036-K: GLP1-118.218 (19 September 2025)

[SigmaAldrich.com/support](https://www.sigmaaldrich.com/support) • [SigmaAldrich.com](https://www.sigmaaldrich.com)

© 2025 Merck KGaA, Darmstadt, Germany and/or its affiliates. All Rights Reserved. The vibrant M logo and MILLIPLEX are trademarks of Merck KGaA, Darmstadt Germany or its affiliates.