

C8721 Cetrinix Supplement

An antibiotic supplement recommended for the selective isolation of the *Pseudomonas* species.

Composition

Cetrimide	100.00 mg
Nalidixic Acid	7.5 mg

Store at 4°C.

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of 0.2N (930-65) sodium hydroxide. Mix well and aseptically add it to 500 ml of sterile, molten (P2102) *Pseudomonas* Agar Base. Mix well and pour into sterile petri plates.

Precautions and Disclaimer

This product is for R&D use only, not for drug, household, or other uses. Please consult the Material Safety Data Sheet for information regarding hazards and safe handling practices.

