3050 Spruce Street, St. Louis, MO 63103 USA Tel: (800) 521-8956 (314) 771-5765 Fax: (800) 325-5052 (314) 771-5757 email: techservice@sial.com sigma-aldrich.com

# **Product Information**

Anti-LC3C antibody produced in rabbit

Affinity isolated antibody

Product Number SAB4200822

### **Product Description**

Anti-LC3C antibody is developed in rabbit using a synthetic peptide corresponding to the N-terminal region of human LC3C (GeneID: 440738), conjugated to KLH as immunogen. The antibody is affinity-purified using the immunizing peptide immobilized on agarose.

Anti-LC3C antibody specifically recognizes human LC3C-I and/or LC3C-II and does not cross react with human LC3A or LC3B. The antibody may be used in various immunochemical techniques including immunoblotting (~16/18 kDa) and immunofluorescence. Detection of the LC3C band by Immunoblotting is specifically inhibited by the immunogen.

Macroautophagy, usually referred to as an autophagy, is a major pathway for bulk degradation of cytoplasmic constituents and organelles. In this process, portions of the cytoplasm are sequestered into double membrane vesicles, called autophagosomes, and subsequently delivered to the lysosome for degradation and recycling.<sup>1,2</sup> Although autophagy is a constitutive cellular event, it is enhanced under certain conditions such as starvation, hormonal stimulation, and drug treatments.<sup>3</sup> Autophagy plays an essential role in cellular differentiation, cell death, and aging and is required for normal turnover of cellular components during starvation. A defective autophagy process may contribute to certain human diseases such as cancer, neurodegenerative diseases, muscular disorders, and pathogen infections.<sup>4,5</sup> Autophagy pathway is evolutionary conserved in all eukaryotic cells.<sup>1</sup>

At least 16 genes encoding for autopahgy (ATG) related proteins that are required for autophagosome formation were identified in yeast by genetic screens. For many of these genes, related homologs have been identified in mammals.<sup>6</sup> The Microtubule-associated protein light chain 3 (MAP1LC3), the mammalian homolog of yeast Atg8/Apg8/Aut7, is essential for the formation of autophagosomes.<sup>7</sup> In addition, three human orthologs of yeast Atg8, named LC3A (MAP1LC3A), LC3B (MAP1LC3B), and LC3C (MAP1LC3C), were identified. These three human isoforms show different expression patterns in human tissues.<sup>8</sup>

The LC3 protein exists in two forms, cytoplasmic LC3-I (18 kDa) and autophagosome membrane associated LC3-II (16 kDa). LC3-II is the result of LC3-I cleavage at the C-terminal glycine by cysteine protease Atg4, which is further modified by the mammalian E1- and E2-like enzymes: Atg7 and Atg3.<sup>7</sup> The amount of LC3-II correlates with the extent of autophagosome formation.<sup>3</sup> Thus, LC3C can be used as an autophagosomal and mitophagy marker.<sup>3,7,9-11</sup>

### Reagent

Supplied as a solution in 0.01 M phosphate buffered saline, pH 7.4, containing 15 mM sodium azide as a preservative.

### Precautions and Disclaimer

For R&D use only. Not for drug, household, or other uses. Please consult the Safety Data Sheet for information regarding hazards and safe handling practices.

# Storage/Stability

Store at –20 °C. For continuous use, store at 2–8 °C for up to one month. For extended storage, freeze in working aliquots. Repeated freezing and thawing is not recommended. If slight turbidity occurs upon prolonged storage, clarify the solution by centrifugation before use. Working dilution samples should be discarded if not used within 12 hours.

# **Product Profile**

<u>Immunoblotting</u>: a working dilution of 1:2000-1:4000 is recommended using whole extracts of HEK-293T cells overexpressing human LC3C.

Immunofluorescence: a working dilution of 1:400-1:800 is recommended using human HeLa cells treated with Bafilomycin A1 for autophagosomes accumulation.

<u>Note</u>: In order to obtain best results in various techniques and preparations, it is recommended to determine optimal working dilutions by titration test.

#### References

- 1. Klionsky, D.J., and Emr, S.D., *Science*, **290**, 1717-21 (2000).
- 2. Kuma, A. et al., Nature, 432, 1032-1036 (2004).
- 3. Kabeya, Y. et al., *EMBO J.*, **19**, 5720-5728 (2000).
- Reggiori, F. and Klionsky, D.J., *Eukaryot. Cell*, 1, 11-21 (2002).
- 5. Shintani, T., and Klionsky, D.J., *Science*, **306**, 990-5 (2004).
- 6. Klionsky, D.J. et al., *Dev. Cell*, 5, 539-545 (2003).
- 7. Tanida, I. et al., *J. Biol. Chem.*, **279**, 47704-10 (2004).
- 8. He, H. et al., J. Biol. Chem., 278, 29278-87 (2003).
- 9. Maier, H.J. et al., Autophagy, 9, 496-509 (2013).
- 10. Le Guerroué, F. et al., *Mol. Cell*, **68**, 786-796 (2017).
- 11. Stolz, A. et al., EMBO J., 36, 549-564 (2017).

VS,DR,OKF,MAM 03/19-1