

shelf life information data sheet

Product Family: Clarisolve® Process Scale Pod

shelf life information

Shelf Life: 60 Months

Expiration Date: N/A

Storage Conditions: Ambient temperature

Applicable Catalog Number(s): CS20MS01F1, CS20MS03F1, CS20MS05F1, CS40MS01F1, CS40MS03F1, CS40MS05F1, CS60HX01F1, CS60HX03F1, CS60HX05F1

quality attributes studied

Visual Inspection

Housing Integrity

Housing Robustness

Not all products listed have undergone stability testing to determine the products' shelf lives. The design of the stability schedule is such that only select product(s) within a product family are tested. The design assumes that the products selected are representative of the family as they utilize the same materials of construction, production processes, and packaging components as the products within the family.

The product(s) above does not have a published expiration date or shelf life claim. The shelf life stated above is a guideline based upon available product stability data which is constantly updated as new data is made available. It is not a guarantee or a claim designating the time prior to which a lot of the product will remain within the shelf life specification if stored under the defined conditions, and after which it must not be put into use. It is not an expiration date or a use-by-date. Further, as this is not a product claim, updates to product shelf lives do not require customer notification.

For firms that must establish expiration dates or use-by-dates for this product(s) utilized in a cGMP manufacturing environment, it is recommended that the firms assign a default use-by-date based on their internal procedures/practices.



16-Oct-19

Thomas Brackett
Supplier Quality Engineering Manager
Third Party Finished Goods Quality

The life science business of Merck KGaA, Darmstadt, Germany
operates as MilliporeSigma in the US and Canada.

